



DMO® Spinal wrap



live your life.

live your life

DMO® spinal wrap



The DMO® Spinal wrap is a simple solution to provide core stability improvement without weight or rigidity.

Provides:

- > Improved trunk control
- > Pelvic and core stability
- > Improved posture
- > Increased sensory feedback

Effective for:

- > Pelvic support
- > Pain relief

Donning and doffing instructions

Detach the Velcro section of the stability wrap, to allow the belt to be opened for donning. The belt would then be positioned on the area of the body requiring support, and the Velcro reattached applying some firm stretch to the material, but ensuring the wrap is comfortable.

To remove the wrap, detach the Velcro.

The DMO® Spinal wrap is often used as an assessment tool for creating the DMO® suit for children.

The DMO® Spinal wrap is available in different lengths and widths. Machine washable. Available in black only.

Sizes:

SSB04 Narrow Fitting

Small	740mm x 80mm
Medium	960mm x 80mm
Large	1140mm x 80mm

SSB03 Wide Fitting

Small	740mm x 160mm
Medium	960mm x 160mm
Large	1140mm x 160mm

Available in black

Nylon/Fluffed finish neoprene

Hook and loop fastening

Intended purpose

The intended purpose of the DM Orthotics Spinal Stability Belt is to provide improved core and pelvic stability without the use of rigid support. The product can be used for pain relief and also can be used by a clinician as an assessment tool or therapy aid in treatment sessions. The compression can also assist with sensory feedback.

Indications

Children or adults who have pain associated with poor alignment or fatigue caused by their inability to maintain a corrected pelvic and lower trunk control. Children and adults who exhibit poor pelvic and core stability and alignment.

Contraindications

Skin conditions that might be exacerbated by close contact with synthetic materials for prolonged periods may require further consultation with a healthcare professional who is overseeing the patient's dermatological treatment, but this is rarely a contra-indication.

Intended Users

The person ordering the product may be the user, a carer, parent or a health care professional.

The wearer of the product will typically be a lay person.

The wearer's carer who may be involved in the donning and doffing of the product, the washing of the belt and the reporting of any issues with the suit will typically be layperson.

SSB03: 5054522000524 (UDI)

SSB04: 5054522000531 (UDI)

DM-003789-FM (Issue 4)

Spinal wrap washing instructions

Before Laundering

- > Ensure all fastenings are closed.

During Laundering

- > Hand wash only in lukewarm water with mild soap or soap flakes.
- > Squeeze garment thoroughly in warm water, rub lightly where soiled.
- > Rinse thoroughly in tepid water.

After Laundering

- > Wrap garment in a towel to absorb moisture.
- > Dry garment flat avoiding direct heat or sunlight.



DO NOT

- > Wash your garment in a washing machine.
- > Dry your garment in a warm air dryer or on household radiators.
- > Wash your garment in bleaching liquids or powders.
- > Leave the garment folded to dry.
- > Use a wringer.

For more information please contact: +44 (0) 1209 219205

www.dmorthotics.com admin@dmorthotics.com

© Copyright 2026DM Orthotics Ltd. Patents and design registrations apply.

