



Spinal Stability Wrap



- Designed specifically to provide spinal stability without weight or rigidity.
- Made from a nylon/fluffed finish neoprene, the wrap is made to wrap around the trunk to provide both stability and proprioceptive feedback for neurologically impaired children.
- Made in two widths and small, medium and large sizes to accommodate infants up to junior children.

MEASUREMENTS		
	Circ (cm)	Deep (cm)
NARROW FITTING		
Small	74	8
Medium	96	8
Large	114	8
WIDE FITTING		
Small	74	16
Medium	96	16
Large	114	16
<i>Made to measure option also available</i>		

For more information please contact +44 (0) 1209 219205
www.dmorthotics.com